

VG Medical

VG Medical is a leading manufacturer for OR and ICU with various product lines of surgical light, operating table, medical supply unit and digital OR system, providing safe, effective and cost-saving solutions for modern hospitals.

Technical Specification

NO.	Instruction	Paramete
1	Length of table top(with head)	2080mm±100mm
2	Width of table top(with side rail)	570mm±50mm
3	Lowest height	700mm±50mm
4	Table top lifting range	≥300mm
5	Trendlenburg/Revers	30°±2°
6	Tilt	20°±2°
7	Kidney elevator	≥110mm
8	Backboard	Upwards ≥75° Downwards 20° ±2°
9	Head board	Upwards ≥45° Downwards ≥50° Detachable
10	Leg board	Downloads ≥90° Outwards≥90° Detachable
11	Load capacity	≥140kg

VG-BCHU-0008.A



COMFORT 200

Electric Operation Table



VG Medical Technology Co.,Ltd.

Tel: +86-510-88816999

Fax: +86-510-88806999

Add: No. 117 GaoYun Road, Wuxi, China

Email: sales@vg-med.com

Web: www.vg-med.com

* Promotional material. Specifications maybe changed without notification.Please contact us for updated information



QUALITY AND STANDARD

- Comply with IEC 60601-2-46, IEC 60601-1, IEC 60601-1-2 ;
- ISO 13485 certified by TÜV;
- MDD 93/42 ;EEC.



PRODUCT FEATURES



- Imported electric motor(actuator). Reliable and easy maintenance. No oil leakage and low noise.
- Fully electric controlled with position sensor. Easy operation with "Reset" ,"Flex" and "Re-Flex" button.
- 125mm low noise medical class castor and central break system. Easy to move and lock firmly.
- SS 304 frame, cover and side rail. Easy cleaning and anti-bacteria.
- Integrated kidney elevator
- Build-in back up battery for supporting more than 100 operations.
- Multi-layer memory foam pad with PU cover. High-frequency welding. Water proof, stain resistant, anti static and easy cleaning.
- High strength image-compatible table board.

SURGICAL POSITION



General Surgery



Cardiac Surgery with CPB



Kidney Surgery



Laparoscope Surgery



Abdominal Surgery



Spinal Surgery - Knee position



Spinal Surgery - Prone position



OB/GYN Surgery