

VG Medical

VG Medical is a leading manufacturer of OR and ICU with various product lines of surgical light, operating table, medical supply unit and digital OR system, providing safe, effective and cost-saving solutions for modern hospitals.



Model and Configuration

Model	COMFORT 300-519	COMFORT 300-519D
Table Length	2050mm	2050mm
Table Width	520mm	520mm
Lowest Height	700±20mm	550±20mm
Height Adjustment	350mm ■	350mm ■
Trendelenburg	30° ■	30° ■
Rev-Trend	30° ■	30° ■
Tilt (Left)	20° ■	20° ■
Tilt (Right)	20° ■	20° ■
Slide	320mm ■	320mm ■
Backboard Upwards	75° ■	75° ■
Backboard Downwards	35° ■	35° ■
Leg board Downwards	90° ▲	90° ▲
Leg board Outwards	90° ▲	90° ▲
Head board Upwards	46° ▲	46° ▲
Kidney Elevator	120mm ▲	120mm ▲
Head Board	Dismountable	Dismountable
Leg Board	Dismountable	Dismountable

COMFORT 300 Load capacity is 280KGS

■ Electrical ▲ Manual



VG Medical Technology Co.,Ltd

Tel: +86-510-88816999

Fax: +86-510-88806999

Add: No.35-201,Changjiang South Road,Wuxi,Jiangsu,China

Email: shanghai@vg-med.com

Web: www.vg-med.com

*Promotional material.Specifications maybe changed without notification.Please contact us for updated information

VG-BCHU-0016.A



COMFORT 300

Electric Operation Table



QUALITY AND STANDARD

- Comply with IEC 60601-2-42, IEC 60601-1, IEC 60601-1-2 ;
- ISO 13485 certified
- MDD 93/42 ;



PRODUCT FEATURES

Germany made electro-hydraulic system, with strong lifting ability and high reliability, free maintenance for three years, dual controlling electric system, intelligent control, recovering with one button, Electric brake, easy to lock

The table top can moved more than 320mm by power-driven, table board is full radiolucent, comply with YY/T0737-2009 stander, it is easy for accessing of C-arm and X-ray, no shooting blind spots

SUS#304 stainless, anti-erosive, easy clean, rubber plastic molded dust cover

Multi-layer memory foam pad, high frequency weld, quick connection device, improve the safety and accuracy

Intelligent battery system with a large capacity battery, it is suitable for at least 100 surgical operations



Germany made electro-hydraulic system



Emergency controlling panel



Integrated kidney elevator



Stainless body and memory foam pad

SURGICAL POSITION

Orthopedic/Trauma surgery



Supine Position (imagery space is 1200mm)



Beach Chair Position



Prone Position (spine surgery)



Supine Position (Orthopedic)

Urology/Gynecology

COMFORT 300 is suitable for different surgical position with the high flexibility



Urology/Gynecology Position



Calculus Position



Lateral Position



Supine Position

General surgery

COMFORT 300 is suitable for surgical and examination with the high flexibility



Supine Position (cardiovascular)



Lateral Position (pelvis and spine surgery)



Supine Position (ENT/Neurosurgical—with head support)



Knee Position



Sitting Position (Neurosurgical—with skull clamp)