

CreLed 150M Examination lamp

CreLed 150M Examination Lamp

Light illumination: > 5000~20000

Lighting depth: > 1250mm

Power supply : 100~240VAC, 50/60Hz

Life of bulb : > 50000 hours

Height of Lamp: 1800mm

Castors: 5castors



Product Attributes

Model No	CreLed 150M
Brand	Lewin
Illumination(LX)	5000~20000
Color Temperature	4000+250
CRI(Ra)	≥90
Dimming Control	4 shift scales
Diameter of Light Head	7cm
Extension of Arm	80cm adjustment
Highest Height of Lamp	1930mm
Power Supply	85-265V AC, 50/60 Hz
Power(W)	5W
Bulb service life more than	50,000 hours
Color of light body	black
Bulb number	1pc

