

CreLed 3400/3300 LED Operation Lamp

German ONDAL arms, German OSRAM bulbs

USA Core Plate, France Lens

680mm/480mm lamp head

lighting intensity of 60000-160000lux

Most of parts imported

excellent focusing and cold light effect

Color temp adjustment with 6 steps from
3500k-5000k

Four modes of background, surface, deep cavity and
normal mode



Product Attributes

Type	CreLed 3400	CreLed 3300
Lamp Head Diameter(mm)	700	600
Quantity of LED Bulbs	64	45
Illumination (Lux)	80,000-160,000lux	60,000-140,000lux
Color Rendering Index (Ra)	≥96	≥96
Color Temperature (K)	3500-5000	3500-5000
Lighting Depth (mm)	≥1200	≥1200
Spot Diameter (mm)	150-350	150-350
Dimming Range	12 levels	12 levels
Operation Mode(Background, Surface, Deep Cavity and Normal Mode)	Yes	Yes
Bulb Service Life (hours)	≥50,000	≥50,000
Camera System	Optional	Optional

