

CreLed 3400/3400 LED Operation Lamp



Product Attributes

Model No	CreLed 3400/3400
Brand	Lewin
Lighting Depth(Mm)	≥1200
Dimming Range	25%-100%,12 Levels
Color Rendering Index (Ra)	≥95
Operation Mode	Back Ground,Surface,Deep Cavity,Normal
Spot Diameter (Mm)	150-350
Lamp Head Diameter(Mm)	700/700
Quantity Of LED Bulbs	64+64
Color Temperature (K)	3500,3750,4050,4350,4650,5000
Bulb Service Life(Hours)	≥80,000
Illumination (Lux)	80,000-160,000

