

CreLed 310M Examination lamp

Operating Light ยี่ห้อ Lewin รุ่น Creled 310M

โคมไฟผ่าตัดเล็กขนาดไม่น้อยกว่า 60,000 ลักซ์

- ปรับความสว่างได้
- ปรับอุณหภูมิแสงได้
- มี Battery

Examination Lamp

Light illumination: $> 60,000\text{lux}$

Lighting depth: $> 1250\text{mm}$

Power supply : $100 \sim 240\text{VAC}, 50/60\text{Hz}$

Life of bulb : $> 50000\text{ hours}$

Height of Lamp: 1-2.8m

Castors: 5 breakable castors

Battery is optional



Product Attributes

Model No	CreLed 310M
Brand	Lewin
light source	imported led bulb
Number of led bulb	14pcs Led bulbs
illumination	$\geq 60,000\text{ lux}$
Color temperaturer	4000,4500,5000
Spot diameter	120-300mm
Depth of lighting	1200mm
Power of LED bulb	1w
Average service life	50000 hours
Certification	CE,ISO9001,ISO13485,FDA
Application	Surgical examination
MOQ	1unit
OEM	Available

THAI MEDICALS TRADING CO.,LTD

67/17 Moo.1 Mahasawat, Phutthamonthon Nakhonpathom 73170 Thailand
www.tmt-med.com,www.thaimedicaltrading.co.th +66 81-3451047 +66 98-6329351



SERVICE

081-345-1047
098-632-9351