

CreLed 3300M Mobile Operation Shadowless Lamp

Mobile Operation Shadowless Lamp

500mm lamp head

lifting height around 1.4m-2.5m

Battery last more than 4hours

lighting intensity of 70000-145000lux

more than 45 LED lights

excellent focusing and cold light effect

with or without battery optional

Four modes of background, surface ,deep cavity and normal mode



Product Attributes

Model No	CreLed 3300M
Brand	Lewin
Lamp Head Diameter	500 mm
Type	Operation Illuminating Lamps
Bulb Number	45 LED Lights
Power Supply	110~130/220~240V, 50/60Hz
Bulb Service Life	≥50,000 Hours
Four Modes	Background, Surface ,Deep Cavity and Normal Mode
Battery Capacity	4~6 Hours
lifting height around	1.4m~2.5m
Illumination	70000-140000 lux

