

This surgical table can be used for different surgery, multifunction use.

For urinary surgery, thoracic surgery, kidney surgery, orthopaedics and other surgery need to use C-arm, X-ray machine, For microsurgical operations, neurosurgery, ophthalmology and otorhinolaryngology, facial/plastic surgery, ENT (ear, throat, eye) surgery spine and kidney surgery, thoracic and renal surgery



It is multifunctional electric operating table with 7 electric adjust functions, electric brake, electric longitudinal slide, electric kidney, electric height adjustment, electric trendelenburg, electric tilt, electric back plate adjustment



- Tabletop Size: 2070mm*525mm
- Height Adjustment(without mattress):680mm-980mm
- Trendelenburg&Reverse Trendelenburg: 30°/ 30°
- Left/Right Tilt: 21°/ 21°
- Back Plate(Up/Down): 82/32°
- Leg Plate adjustment (Up/Down): 50°/93°
- LegPlate(Horizontal adjustment): 90°
- Electric built-in kidney bridge :126mm
- Longitudinal slide: 310mm
- Load capacity: 270kg

- Double control system,Column panel controller and wire hand controller
- One key reset
- Head and leg plate interchange
- High-strength carbon steel table frame,
- 304 stainless steel base and column
- "T" shaped base more stable,convenient for doctor
- Mattress material : memory sponge foam 70mm thickness, can be X-ray perspective
- Column emergency : can last 2 mins to adjust position when AC power off
- Position memory: 3 sets, can be saved and deleted
- C-arm & X-ray compatible
- Can compatible with orthopedic traction frame







Electric Lock



Back Plate Up 82°



Back Plate Down 32°



Height Up 680mm



Height Down 80mm



Leg Plate Up 50°



Leg Plate Down 93°



Left Tilt 21°



Right Tilt 21°



Longitudinal slide 310mm



Trendelenburg 30°



Reverse Trendelenburg 30°



Electric bridge 126mm



Electric bridge 0 mm

Standard Accessories



Table base cover

The base material for the table is high thickness 304 stainless steel, seamlesed welding by robot, with modern appearance and more stable.



Mattress

Memory and high density sponge mattress , waterproof and seamless , avoid liquid in . the connection between table and mattress is buckle joint , easy to disassemble and install , easy to clean.



Head plate

Head plate is with imported gas spring for easy control.



Kidney bridge

Liftable electric kidney bridge.

Horizontal sliding

Extra horizontal sliding , distance 310mm
Memory sponge mattress ,
C-arm and X-ray compatible



Dual control system

All functions are controlled by both hand controller and control panel on table base column.
with control panel on table base column in case of hand controller failure

Leg plate

Leg plate with gas spring for easy control ,
disassembled type leg plate





Arm plate

Anesthesia screen



Leg holder

The special structure which can prevent the shaking during the operation

Shoulder support

The distance and angle
can be adjusted



Waist support

The distance and angle
can be adjusted