



Shenzhen Haiye Medical Technology Co., Ltd.



HYHJ-KC Reusable Video Laryngoscope

Uses and Application

It is used to clinically provoke the epiglottis to expose the glottis, and guide the medical staff to perform tracheal intubation and mechanical ventilation under visible conditions.

Applicable to the anesthesiology department, ICU, pediatrics, emergency department, ambulance, etc..

➤ Product Features and Configuration

The monitor of this product can be adapted 6 sizes camera system (blades): Miller0, Miller1, Mac2, Mac3, Mac4, D-blade

This product has the functions of taking pictures, recording, freezing, data storage and export, picture and video playback, HDMI, Type-C interface, etc., and can also be connected to an external screen to display on the same screen.



➤ Basic Parameters

Screen	3.5" color display
Camera	Pixel:3M; Field depth:6-100mm; Field angle:70° ±9° ; Resolution:7.13 lp / mm
Battery	Working time:6h; Charging time:4h; Capacity:3400mAh
Function	HDMI, Type-C, support photo, video, freeze frame, storage and data export
Storage	Environment temperature -40°C~70°C
Humidity	≤ 93% rh



Tel: 0086-029-68069218
Phone No.: 0086-18629338051
Email: HYYL2013@haiyemedical.com
Website: www.hyendoscope.com