

Hollow CreLed 5500/5500 LED Operation Lamp

German ONDAL arms, German OSRAM bulbs

530mm lamp head

lighting intensity of 60000-120000LX

more than 21 LED bulbs:

excellent focusing and cold light effect

Electric Spot adjustment function

Shadow control system

Color temp. 3500k-5000k

Four modes of background, surface, deep cavity and normal mode



Product Attributes

Type	Hollow CreLed 5500/5500
Lamp Head Diameter(mm)	500
Quantity of LED Bulbs	21
Illumination (Lux)	60,000-120,000
Color Rendering Index (Ra)	≥96
Color Temperature (K)	4000±500
Lighting Depth (mm)	≥1200
Spot Diameter (mm)	150-350
Dimming Range	12 levels
Operation Mode(Background, Surface, Deep Cavity and Normal Mode)	Yes
Bulb Service Life (hours)	≥50,000

