









SonoScape



S11 Plus

Surpass Your Expectations

SonoScape

SonoScape Medical Corp.
2F, 12th Building, Shenzhen Software Park Phase II, Keji Middle
2nd Road, Shenzhen, China
Tel: +86-755-26722890 Fax: +86-755-26722850
Email: market@sonoscape.net www.sonoscape.com



U-S11 Plus20190527

S11 Plus


- another powerful expression of SonoScape's trolley color doppler system, brings you an upscale experience in its class. Quality imaging embodies enhanced exam consistency, helping you diagnose assuredly through versatile applications and attention to details. With improved workflow and intelligent features, S11 Plus is the system you can rely on with ease.

ENHANCED IMAGING

- PIH - Pulse Inversion Harmonic Imaging




- SCI - Spatial Compound Imaging



- Multi-beam Forming Technology



- μ -Scan – Speckle Reduction & Edge Enhancement



μ -Scan off



μ -Scan on

STRENGTHENED APPLICATION

- Versatile Transducer Family
- Professional Measurement Package
- 3D/4D Imaging
- Panoramic Imaging
- Vis-Needle
- Anatomical M Mode

SIMPLIFIED WORKFLOW

- Auto IMT
- Auto EF
- Auto NT
- Auto Trace
- Auto Optimization

INGENIOUS ERGONOMIC DESIGN

- High-definition medical LED monitor with articulating arm
- Four interoperated probe sockets for wide application
- Built-in battery for 90 min continuous scanning
- Simplified user-interface for efficient operation