



Smart◎Large◎Rapid



2024-Catalog



ZHEJIANG GETIDY MEDICAL INSTRUMENT CO., LTD.
ADD.: DUNCUN, DONGNAN INDUSTRIAL ZONE, SHUXI STREET, WUYI, JINHUA,
ZHEJIANG 321200, CHINA.
Tel: +86-579-87712106 Web: www.getidy.com www.getidydental.com
212725035 +86 15705893636
E-Mail: john@getidy.com getidy@getidy.com Mobile Phone: +86 15705893636

2024.11.30



PROFESSIONAL MANUFACTURING
STEAM STERILIZER

Steam sterilizer >>> CLASS B


Model:YS-29 YS-45



European Class B and S, in accordance with EN 13060-2014 standard, three times(CLASS B) and one time(CLASS S) pre-vacuum sterilization.Vacuum index: -85Kpa.Residual humidity <0.2%.Being able to sterilize instruments those are packed, unpacked,solid,Hollow A,Hollow B and Hollow B pipeline.High efficiency,all wave sterilizing.

- ※ 800*480 Colour touch Screen
- ※ Automatic door structure
- ※ Water Monitoring
- ※ User-defined Function
- ※ New Rapid Steam Generator
- ※ Heraeus Temperature Sensor(dual)
- ※ USB Records
- ※ Optional:External Printer
- ※ Resolution is 1kPa and 0.1℃
- ※ Delay-Start/Sleep-Mode
- ※ Pipeline Cleanning
- ※ Helix and BD Test
- ※ MXP 5700 Pressure Sensor
- ※ Build-in Printer

Model	Voltage(V/Hz)	Power(VA)	Vol(L)	Chamber(mm)	Product Size(mm)	Pack Size(mm)	G.W/N.W(kg)	Space(cm)
YS-29	230/50/60	5500	29	Φ319×420	688×565×560	835×705×810	126.5/109.0	70
YS-45	230/50/60	5500	45	Φ319×620	838×565×560	985×705×810	140.0/121.5	80



>Wide-Color Touch Screen



>Download Sterilization Records by U disk



>Heating Water Pump



>Double condensers and four fans



Model:YS-18 YS-23



European Class B and S, in accordance with EN 13060-2014 standard, three times(CLASS B) and one time(CLASS S) pre-vacuum sterilization.Vacuum index: -85Kpa.Residual humidity <0.2%.Being able to sterilize instruments those are packed, unpacked,solid,Hollow A,Hollow B and Hollow B pipeline.High efficiency,all wave sterilizing.

- ※ 480*272 Colour Touch Screen
- ※ Delay-Start
- ※ Sleep-Mode
- ※ Pipeline Cleaning
- ※ Build-in Printer
- ※ Optional: Water Monitoring
- ※ MXP 5700 Pressure Sensor
- ※ Heraeus Temperature Sensor(dual)
- ※ Air Pressure Control Lock
- ※ Multi Languages
- ※ Resolution is 1kPa and 0.1℃
- ※ External Printer

Model	Voltage(V/Hz)	Power(VA)	Vol(L)	Chamber(mm)	Product Size(mm)	Pack Size(mm)	G.W/N.W(kg)	Space(cm)
YS-18	230/50/60	1600	18	Φ250×355	565×460×480	635×535×535	55.0/64.0	60
YS-23	230/50/60	1600	23	Φ250×450	565×460×480	635×535×535	56.5/65.5	60

Model:YG-18 YG-23



European Class B and S, in accordance with EN 13060-2014 standard, three times(CLASS B) and one time(CLASS S) pre-vacuum sterilization.Vacuum index: -85Kpa.Residual humidity <0.2%.Being able to sterilize instruments those are packed, unpacked,solid,Hollow A,Hollow B and Hollow B pipeline.High efficiency,all wave sterilizing.

- ※ White LED
- ※ MXP 5700 Pressure Sensor
- ※ Heraeus Temperature Sensor
- ※ New Rapid Steam Generator
- ※ USB Records
- ※ Optional:External Printer
- ※ Ulka Water Pump
- ※ Pipeline Cleaning
- ※ User-defined Function
- ※ Helix and BD Test
- ※ Resolution is 1kPa and 0.1℃

Model	Voltage(V/Hz)	Power(VA)	Vol(L)	Chamber(mm)	Product Size(mm)	Pack Size(mm)	G.W/N.W(kg)	Space(cm)
YG-18	230/50/60	1600	18	Φ250×355	575×460×445	635×535×490	60.0/50.5	60
YG-23	230/50/60	1600	23	Φ250×450	575×460×445	635×535×490	61.0/52.0	60

Model:YY-18 YY-23



European Class B and S, in accordance with EN 13060-2014 standard, three times(CLASS B) and one time(CLASS S) pre-vacuum sterilization.Vacuum index: -85Kpa.Residual humidity <0.2%.Being able to sterilize instruments those are packed, unpacked,solid,Hollow A,Hollow B and Hollow B pipeline.High efficiency,all wave sterilizing.

- ※ White LED
- ※ MXP 5700 Pressure Sensor
- ※ Heraeus Temperature Sensor
- ※ New Rapid Steam Generator
- ※ USB Records
- ※ Optional:External Printer
- ※ Ulka Water Pump
- ※ Pipeline Cleaning
- ※ User-defined Function
- ※ Helix and BD Test
- ※ Resolution is 0.1bar and 1℃

Model	Voltage(V/Hz)	Power(VA)	Vol(L)	Chamber(mm)	Product Size(mm)	Pack Size(mm)	G.W/N.W(kg)	Space(cm)
YY-18	230/50/60	1600	18	Φ250×355	575×460×445	635×535×490	60.0/50.5	60
YY-23	230/50/60	1600	23	Φ250×450	575×460×445	635×535×490	61.0/52.0	60