

Technical Specifications

Model	ZW-500D
Lamp Head Diameter	500mm
Light Source	Imported Osram
LED Quantity	48pcs
Power	65VA
Work Voltage	AC 100-240V
Frequency	50Hz/60Hz
CRI (Ra)	≥85
CCT (Color Temp)	4000K-5000K (adjustable)
Brightness	40000lux -1200000lux
Spot Diameter	120-280mm
LED Life	50,000hours
Illumination depth	≥800mm
Temp. on Surgery Head	≤1 °C
Standard Install Height	2800-3200mm



Wall Mounted or Ceiling Type for Options