


Technical Specifications

Model	ZW-700/500D	
Lamp Head Dia.	700mm	500mm
Light Source	Imported Osram	Imported Osram
LED Quantity(pcs)	80	48
Power	70VA	65VA
Work Voltage	AC 100-240V	
Frequency	50Hz/60Hz	
CRI (Ra)	≥85	
CCT (Color Temp)	4000K-5000K (adjustable)	
Brightness (0-100stages adjust)	40000lux - 1600000lux	40000lux- 120000lux
Spot Diameter	120-260mm	120-280mm
LED Lifespan	50,000hours	50,000hours
Illumination depth	≥600mm	≥800mm
Temp. on Surgery Head	≤1 °C	≤1 °C
Install Height	2800-3200mm	2800-3200mm
Others	CE, ISO13485, Free Sales Cert.	



Hunan Mplent Optoelectronics Technology Co.,Ltd.

Email: jena@mplent.net Cell: +86 18175116933