

iMAC10

*ECG Acquisition Unit for
Telemedicine*

zoncare

Wuhan Zoncare Bio-medical Electronics Co.,Ltd.

Add: Zoncare Building, No. 380, High-tech 2nd Road, Eastlake High-tech
Development Zone, Wuhan, Hubei 430206 P. R. China.
Tel: +86-27-8777 0203 Fax: +86-27-8777 0581


Email: info@zoncare.com <http://www.zoncare.com>

2021 Wuhan Zoncare Bio-Medical Electronics Co.,Ltd. All rights reserved. Specifications subject to changes
without prior notice.

zoncare


Lead connection diagram

- Support lead-off sound and light alarm
- Support battery reminder alarm
- Flexible connectivity
- Bluetooth connection
- Idea for Telemedicine Application
- 12 Simultaneous Lead Acquisition and Display
- The sampling rate is 32000hz to make the waveform restoration more accurate
- Able to detect pacemaker signals to assist doctor in clinical diagnosis



Ingenious Design

External USB magnetic connection charging, safe and convenient




Zoncare provides a proprietary wireless sampling module ,bringing a new ECG experience with familiar resting ECG operations on both smart phone, tablet and PC.

Data Transmission Resting ECG App(for Android)

- Real Time Sampling
- Automatic Diagnosis
- Touch screen can provide a better user experience

Data Transmission Workstation(for PC)

- Printing
- Waveform Analysis
- Average Template
- Automatic Diagnosis
- Measure
- Reanalysis



Optional Trolley

Filter

Armed with the comprehensive filters, users could acquire superb ECG waveform with convenient ways .

ECG Sampling

Both pre-sampling and real-time sampling modes are provided to meet different situation, no matter it is an emergent case, home visiting, or general ward rounds.